

A Guide to Matching Filter Tips



Brand	Specifications		Filter tips										
			10ul Short	10ul Long	20ul	100ul	200ul	200ul Long	300ul	1000ul	1000ul long	1250ul	
Thermo	Single Channel	F3 1-10ul	o	o									
		F3 2-20ul			o	o	o	o	o				
		F3 5-50ul				o							
		F3 10-100ul				o							
		F3 20-200ul				o	o	o	o				
		F3 100-1000ul								o	o	o	
		F2 0,5-15ul											
		F2 1-10ul	o	o									
		F2 0,2-2ul	o	o									
		F2 10-100ul				o							
	F2 100-1000ul								o	o	o		
	Multi Channel	F2 1-10ul	o	o									
		F1 1-10ul	o	o									
F3 30-300ul													
Sartorius	Single Channel	proline 2-20ul			o	o							
		proline 10-100ul				o							
		proline 20-200ul				o	o	o	o				
		proline 100-1000ul								o	o	o	
	Multi channel	0,5-10ul	o	o									
Biohit	Multi Channel	20-200ul					o	o	o				
Brand	Single Channel	0,1-2,5ul	o	o									
		1-10ul	o	o									
		2-20ul			o	o							
		10-100ul											
		20-200ul				o	o	o	o				
	100-1000ul								o	o	o		
	Multi Channel	30-300ul							o				
Gilson	Single Channel	PN 2-20ul			o	o							
		PN 10-100ul				o							
		PN 0,2-2ul	o	o									
		PN 20-200ul					o	o	o				
		PN 100-1000ul								o	o	o	
Eppendorf	Single Channel	0,1-2,5ul	o	o									
		0,5-5ul	o	o									
		1-10ul	o	o									
		2-20ul			o	o							
		10-100ul				o							
		20-200ul					o	o	o				
	100-1000ul								o	o	o		
	Multi Channel	0,5-10ul	o	o									
		5-100ul				o							
		50-100ul				o							
		20-200ul					o	o	o				
30-300ul								o					
50-1200ul													
Dragonlab	Single Channel	0,1-2,5ul	o	o									
		0,5-10ul	o	o									
		1-10ul	o	o									
		2-20ul			o	o							
		10-100ul				o							
		20-200ul					o	o	o				
	100-1000ul								o	o	o		
	Multi Channel	10ul	o	o									
300ul							o						

A Guide to Matching Filter Tips



Brand	Specifications		Filter tips									
			10ul Short	10ul Long	20ul	100ul	200ul	200ul Long	300ul	1000ul	1000ul long	1250ul
Rainin	Single Channel	L10	○	○								
		L-20XSL, 2-20ul	○	○								
		L-200PL+, 20-200ul					○	○				
		L-1000PL+, 100-1000ul								○	○	○
		SL-200XSL, 20-200ul					○	○				
		SL2, 0.1-2ul	○	○								
		SL200, 20-200ul					○	○				
		SL1000, 100-1000ul								○	○	○
Terch	Single Channel	0.1-2ul	○	○								
		0.5-10ul	○	○								
		2-20ul			○	○						
		5-50ul				○						
		10-100ul				○						
		20-200ul					○	○	○			
		100-1000ul								○	○	○

Our pipette tip compatibility chart is constantly updated with the latest information, but if you ever find that a product you're looking for is missing or you have any questions, our helpful support team is always available to assist you. Don't hesitate to reach out to our helpdesk for any additional information you may need.