

A Guide to Matching Pipette Tips



Brand	Specifications		Pipette tips								
			10ul Short	10ul Long	200ul	200ul Long	300ul	1000ul	1000ul long	1250ul	
Thermo	Single Channel	F3 1-10ul	○	○							
		F3 2-20ul			○	○					
		F3 5-50ul			○	○					
		F3 10-100ul			○	○					
		F3 20-200ul			○	○					
		F3 100-1000ul							○	○	○
		F2 0,5-15ul	○	○							
		F2 1-10ul	○	○							
		F2 0,2-2ul	○	○							
		F2 10-100ul			○	○					
	F2 100-1000ul							○	○	○	
	Multi Channel	F2 1-10ul	○	○							
		F1 1-10ul	○	○							
F3 30-300ul											
Sartorius	Single Channel	proline 2-20ul			○	○					
		proline 10-100ul			○	○					
		proline 20-200ul			○	○	○				
	Multi Channel	proline 100-1000ul						○	○	○	
Biohit	Multi Channel	0,5-10ul	○	○							
		20-200ul			○	○	○				
Brand	Single Channel	0,1-2,5ul	○	○							
		1-10ul	○	○							
		2-20ul			○	○					
		10-100ul			○	○					
		20-200ul			○	○	○				
	Multi Channel	100-1000ul						○	○	○	
	Gilson	Single Channel	PN 2-20ul			○	○				
PN 10-100ul					○	○					
PN 0,2-2ul			○	○							
PN 20-200ul					○	○					
PN 100-1000ul									○	○	○
Eppendorf	Single Channel	0,1-2,5ul	○	○							
		0,5-5ul	○	○							
		1-10ul	○	○							
		2-20ul			○	○					
		10-100ul			○	○					
		20-200ul			○	○	○				
	Multi Channel	100-1000ul						○	○	○	
	Multi Channel	0,5-10ul	○	○							
		5-100ul			○	○					
		50-100ul			○	○					
		20-200ul			○	○	○				
30-300ul				○	○	○					
Dragonlab	Single Channel	0,1-2,5ul	○	○							
		0,5-10ul	○	○							
		1-10ul	○	○							
		2-20ul			○	○					
		10-100ul			○	○					
		20-200ul			○	○	○				
	Multi Channel	100-1000ul						○	○	○	
Multi Channel	10ul	○	○								
	300ul						○				

A Guide to Matching Pipette Tips



Brand	Specifications		Pipette tips							
			10ul Short	10ul Long	200ul	200ul Long	300ul	1000ul	1000ul long	1250ul
Rainin	Single Channel	L10	○	○						
		L-20XSL, 2-20ul	○	○						
		L-200PL+, 20-200ul			○	○				
		L-1000PL+, 100-1000ul						○	○	○
		SL-200XSL, 20-200ul			○	○				
		SL2, 0.1-2ul	○	○						
		SL200, 20-200ul			○	○				
		SL1000, 100-1000ul						○	○	○
Terch	Single Channel	0.1-2ul	○	○						
		0.5-10ul	○	○						
		2-20ul			○	○				
		5-50ul			○	○				
		10-100ul			○	○				
		20-200ul			○	○	○			
		100-1000ul						○	○	○
Axygen	Single Channel	AP-2-P 0.1-2ul	○							
		AP-10-P 0.5-10ul	○							
		AP-20-P 2-20ul			○					
		AP-50-P 5-50ul			○					
		AP-100-P 10-100ul			○					
		AP-1000-P 100-1000ul						○		
	Multi Channel	AP-8-10-P 0.5-10ul	○							
		AP-12-10-P 0.5-10ul	○							
		AP-8-50-P 5-50ul			○					
		AP-12-50-P 5-50ul			○					
		AP-8-200-P 20-200ul			○		○			
		AP-12-200-P 20-200ul			○		○			
		AP-8-300-P 30-300ul					○			
		AP-12-300-P 30-300ul					○			

Our pipette tip compatibility chart is constantly updated with the latest information, but if you ever find that a product you're looking for is missing or you have any questions, our helpful support team is always available to assist you. Don't hesitate to reach out to our helpdesk for any additional information you may need.